MD

STERILE EO

Sunmatic The Safety Company AB

Sörredsvägen 113 | 418 78 Göteborg | Sweden
+46 31 799 16 00 | info@housegard.com
www.housegard.com

Contains CE-approved
medical products.

RFX+CARE®

RFX+CARE Manufacturing Co., Ltd.
No. 599 Wuxie Road Yuecheng District Shaoxing,
Zhejiang, China 312000

EU REP
RFX+CARE International A/S
Gammel Strandvej 22B, 2990 Nivå, Denmark

STERILE WOUND DRESSING LARGE x1 17 x 17 cm

Intended use: Bandage with compress for protecting a clean and dry wound.

Instructions of use: Clean and dry the wound and the skin around the wound. Apply the compress on the wound and wrap the bandage around. Finish by tying the two ends together.

Warning: Replace the dressing on a daily basis, or more often if necessary. Single use only. Check the integrity of the individual sachet before use. Store in a cool, dry place. Use before the expiry date on the individual sachet. Must not get wet.

ELASTIC GAUZE BANDAGE x1 8 cm x 3 m

Intended use: To hold compresses or other medical dressings in place.

Instruction for use: Wrap the bandage around the area where compress/dressing is placed. Do not apply excessive pressure. Place the end under the lower layers or use a safety pin.

Warning: Never apply the bandage directly onto injured skin. Avoid pressing too hard. In case of discomfort, release the bandage to reduce the pressure. Store in a cool, dry place.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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